

SV Transformed cells for deter.  
mRNA annealing & fragments

1. SV-B3T3 Cl 6A

2. SV-B3T3 Cl 5A

3. UV-15 Cl 5

4. T2

5. SV 3T3 9/8a

6. SV 3T3 Cl 6



Antibody ag 1) infected BSC-1  
2) SV-B-373

Idea: To identify SV40 proteins in productively infected & transformed cells, ppt. the  $^{35}\text{S}$  labeled proteins in these cells & antisera made ag. inf. BSC-1 (early & ara C + late) & SV-Ball 373. Each serum is to be absorbed & using or untr. cell extract. Imm. precipitate is washed & taken up & SDS-ME-EDTA & electrophoresed. Control = using or untr. cell extr. ppted & same antiserum.

Cell extracts for immunization:

1. Normal BSC-1
2. Inf. BSC-1 - ara C = early pr.
3. Inf. BSC-1 - 48 hrs = early + late
4. Normal Ball/373
5. SV transformed Ball/373 - T2
6. " " " " " "

Note: If all genes are expressed in #3 this serum may be okay for testing all other cell extracts!



## Liver RNA polymerase

SHEET NO.

BY

DATE

SUBJECT

30 g. of liver from  
 washed & minced in sol'n A at 0°  
 + 60 ml sol'n A - homogenized in P-E in Teflon pestle.  
 10-15 strokes at ~2000 rev/min.

Filtered through cheesecloth.

Volume made to 150 ml in sol'n A & layered  
 onto sol'n B - 25 ml onto 5 ml B in SW 25 rotor.

Cent at 22,000 rpm for 1 hr x 4. Pour off super & upper tube.  
 Suspend pellets in total of 10 ml of E & sonicate  
 for 6-15 sec. periods. Cent at 80,000 g for 45'.



9/16/69 Cellulose phosphate

SHEET NO.

BY

DATE

SUBJECT

Whatman # 21111 P11  
Medium fibrous powder  
Nominal total capacity 7.4 meq/g

Washed  $\bar{c}$  0.5 N NaOH water, HCl water  
Equil.  $\bar{c}$  0.05 M Tris HCl pH 7.8



9/15/69

## Liver RNA polymerase

SHEET NO.

BY

DATE

SUBJECT

- A. 0.25 M sucrose, 0.05 M TrisCl 7.5,  
0.025 M KCl, 0.005 M  $MgCl_2$  (TKM)  
0.0015 M  $CaCl_2$
- B. 2.3 M sucrose in TKM +  $Ca^{++}$
- C. 2.2 M sucrose, 0.001 M  $MgCl_2$
- D. 0.32 M sucrose 0.001 M  $MgCl_2$ , 0.02 M Tris 7.5
- E. 0.05 M TrisCl 7.8, 30% glycerol, 0.005 M ME

Stocks:

5M sucrose -  
1M Tris 7.5 - 500 ml  
1M Tris 7.8 - 500 ml  
1M  $MgCl_2$  - 100 ml  
1M KCl - 1000 ml

1M = 34.2%

230

2.5

342  
23  
1026  
684  
786.68/2  
1573.38/200

342  
23  
1026  
684  
786.68/2  
1573.38/200

342  
23  
1026  
684  
786.68/2  
1573.38/200



10/28/69

SUBJECT

3T3 }  
SV 3T3 } for KK Takemoto

1 small flask SV 3T3

1 large " 3T3

Med. changed 10/27 MEM = 15% FBS

Next day trypsinized @ 0.5% tryp in PBS - 30'

Centrifuge + wash @ 1 ml med.

Suspend SV 3T3 in 4 ml

3T3 in 8 ml

Take 0.8 of each + 4 ml med to each of 2 flasks

Remainder + 0.6 + 1.6 50% glycerol, resp.

+ put at 4° overnight (10/28-29)



BBL, Cat. #40602, Anti-Rabbit Globulin, Fluorescein labeled, Lot #9061907, has been tested in an indirect staining system employing Salmonella "O" Group D antiserum prepared in rabbits and a Salmonella typhi antigen. Satisfactory results were demonstrated at a conjugate dilution of 1:40.

These results were obtained employing a Zeiss Standard RA binocular Fluorescent Microscope, equipped with a BG-12 exciter filter and a Zeiss 50 eyepiece barrier filter. An HBO-200 mercury burner in a Zeiss housing served as the light source.

It is suggested that each laboratory determine the optimal staining titer under its own standard operating procedures.

BBL, Division of BioQuest  
Cockeysville, Maryland



12/11/69 SV40 DATA f. SVinf CV1

SHEET NO.

BY

DATE

SUBJECT

9-100mm dishes f. Wc - confluent CV1  
 Trypsinizing + suspend cent cells in  
 200 ml MEM ± 10% FBS  
 (Cell count = 64 per 8 small squares  
 = 80/cu mm or  $1.6 \times 10^7$  total or  $8 \times 10^5$ /10 ml)  
 Dispense 10 ml into dishes  
 15 ml into 1 flask

(Orig pl. : had  $\frac{1.6 \times 10^7}{9} = \sim 2 \times 10^6$ /plate)

12/12 Cells nearly confluent

12/13 Wash cells once ± PBS  
 ( $5 \times 10^8$  f.u./ml SP)

3  $\times 10^6$   
 3  $\times 10^7$   
 1.5  $\times 10^8$   
 10:45 AM I - Infect ± SV40 0.15 ml of - 16 plates  
 moi = 30  
 Un - Add 0.15 ml medium - 2 plates

10:45 AM after  $\frac{1}{2}$  hrs at 37° add 10 ml regular MEM ± 10% FBS  
 to group I-A - 10 plates  
 To I-B+C add MEM ± 10% dial. FBS - 10 plates  
 To Un add 10 ml " " " "

12/14 Change all media except I-A

10:45 AM → Add O-thymidine to I-B + Un-B 1:5 = 15 ml  
 24 hrs 100  $\mu$ C/ml 50  $\mu$ C/205 mg (3  $\mu$ C)  
 Add H-thymidine to I-C + Un-C .05/ml  
 1 mC/ml 5 mC/0.675 mg NETO27X CH<sub>3</sub> lab (50  $\mu$ C)



I-A inf - no radioactivity  
 I-B inf +  $^{14}\text{C}$ -TdR  
 I-C inf +  $^3\text{H}$ -TdR  
 Un-B uninf +  $^{14}\text{C}$ -TdR  
 Un-C uninf +  $^3\text{H}$ -TdR

2 hrs.  
 Remove med. + wash cells  $\times 2$   $\bar{c}$  5 ml  
 Tris-saline (per L 0.1 g  $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ ; 0.1 g  $\text{CaCl}_2$ ;  
 8 g NaCl; 0.38 g KCl; 0.1 g  $\text{Na}_2\text{HPO}_4$ ;  
 3 g Tris - adj to pH 7.4)

To each dish add 1 ml 0.6% SDS - 0.001 M EDTA  
 pH 7.5. after 10-20' at rm temp

scrape lysate  $\bar{c}$  rubber policeman + pour  
 into plastic cent. tube - 8 mm diam.

Add 5 M NaCl to  $\rightarrow$  1 M + slowly  
 invert 10 x. store at  $4^\circ$  for  $\geq 8$  hrs.

Cent. at 12500 rpm ( $\sim 17,000g$ )  $\times 30'$  in cold.

Super removed  $\bar{c}$  pipette tip to glass tubes. Keep pellets in refrig.

Count 50  $\mu$ l of each super except A - filter paper  $\bar{c}$   
 Cold thymidine +  $\gamma$ -RNA carrier. Cold TCA wash.  
 Store super in refrig overnight.



12/17

SV40 DNA con't.

SHEET NO.

BY

DATE

SUBJECT

Supers extr. 2X  $\bar{c}$  phenol out  $\bar{c}$  1M Tris HCl pH 8.0  
 then 1X  $\bar{c}$  chloroform - isoamyl alc (24:1).  
 (Sample C lost in centrifuge) -

Supers dial ag 1X SSC  $\bar{c}$  .001M EDTA .01M Tris 7.4  
 Overnight in cold.  
 .010 ml B taken for cts.

### CsCl run

	Tube 1	2	3
	A	A	B
Wt sol'n	2.7584 g	2.9967	3.1928 g
Eth Br	0.3 of 200x/ml	0.3	0.3
CsCl	2.89 g	3.11 g	3.30 g
CsCl sol'n	0.75 ml	0.46 ml	0.21 ml
48.58% w/w			

Dens. 1.572 by weight

Adjusted by adding .02 ml water each tube

Spin in SW50.1 at 45000 rpm x 48 hrs. at 34°F Temp. setting

### From Table

ref index	1.3860	=	1.5522 g/ml	=	48.00% w/w
	1.3891	=	1.5874 g/ml	=	50.00 "
by interp	<u>1.3869</u>	=	1.565	=	48.58 "



1.3869

1.572  
1.566  
— .006

7.10  
.0038  
.027

%  
by wt.  
48.00

1.3860 - 743.7 g/L

1.5522

1.3891 -

792.3

1.5874

50.00

0031

.0352

2.00

0.58

1.3869

100.00  
48.58g  
51.42

Water

2.89  
.234  
— 3.124

3.11  
.234  
— 3.344

.01  
3.30  
.234  
— 3.534

1.572

1.5660

(1.566)

1.5522

.0138

~~Trial - 13.2358~~ 3.5

Beaker 27.0402

.75

48

1.17

5.95

7.12

3/0.7350

27.0402

27.7752

30.9680

4.85  
.234  
— 5.084

(3) # 1 0.7422  
+ 27.0402  
— 27.7824

+ 30.5408

2.7584

+ .3

3.0584

Coll 2.89

2 0.7421

27.0402

27.7823

(7.12) 30.7790

5.95 2.9967

+ .3

3.2967

3.11

(.46)

7.12 6.40

.72

3.1928

+ .3

3.4928

3.30


7.12  
6.79  
— .33  
.21



12/21

SUBJECT

CsCl run - each tube had 2 red bands  
each fluorescent  $\bar{c}$  UV light



I & II collected by drops (A) or by  
past pipette (B).

Eth. br. extr.  $\bar{c}$  isopropanol  $\sim$  equal vol X 3.  
+ DNA sol'n dial ag 2 L  $\frac{1}{10}$  SSC  $10^{-4}$  M EDTA X 2.  
in cold.

Count. 0.25 ml  $^{14}$ C DNA fx.

		Total vol	Total cts
I	416 cpm	$\sim .7$ ml	$\sim 28000$
II	625 "	$\sim 1$ ml	
bet I & II	33 "		



12/20

Cleavage of SV40 DNA by cell extracts

SHEET NO.

SUBJECT

Prep. of extracts (Sambrook + Fritschy J Virol 1969)

BY

DATE

Cells: 1. CV-1 MEM 10% FBS to confluence. 12/20

Fed 24 hrs  $\bar{a}$  washing + freezing

✓ 2. SV 3T3 as above confl.

3. 3T3

✓ 4. MA 196 - Human skin <sup>diplont</sup> f. W. Carter -  $\sim 2/3$  confluent

2/5/70 ✓ 5. L cells - (f. W. Carter) confl (2/3/70) - 5 dishes (100 mm)

Extract Wash cells on dish  $\bar{c}$  ice cold PBS x 2,  
then TED (0.02M TrisCl pH 8, 0.001 EDTA, 0.5mM DTT) x 1.Collect  $\bar{c}$  rubber policeman & store at  $-70^{\circ}$  in  
0.2 ml batches. Before assay thaw & cent at 7000g  
for 10' at  $0^{\circ}$ . Inc supernatantLigase assay conditions:0.1M KCl, 0.04M Tris 7.7, 0.01M  $MgCl_2$ 0.01M ME  $10^{-4}$ M ATP Vol = 0.12 ml

5 of 2 (2/105 prot.) Stopped by add'n EDTA.

1/2



12/29/69 Hemophiles R env. on SV40<sup>14</sup>C-DNA  
for electrophoresis

SHEET NO.

BY

DATE

SUBJECT

Buffer mix 10X

10ml

Tris HCl 7.4 0.10 M

1ml 1M

Mg Cl<sub>2</sub>

0.09 M

.9ml 1M

ME

0.07 M

.05ml 14M

NaCl

0.4 M

.8ml 5M

7.25 water

1

2

Buffer mix

.005

.005

Water → .05

.005

.003

DNA I<sup>14</sup>C

.04

.04

(≈ 800 cpm)

prep B 12/69

Enzyme

0

.002

f. AS 9/12/69 10+12

35° X 30'

Then add 20μl 0.2M EDTA pH 7.4 + .01% SDS

Electrophoresis entire amt. → add'n 5μl 2.5M sucrose = BPB  
(Note: BPB blue in tube 1, yellow in tube 2!)

Buffer Tris - NaAc - EDTA pH 7.8 ± 0.2% SDS

Tube 1 - gel 2

Tube 2 - gel 5

4cm Amp/Tube

6<sup>15</sup> PMto 7<sup>45</sup> PM13<sup>1/2</sup> hrs

5% gels - 10 cm length

~~Results~~ Gels frozen at -70° & sliced in 3 segments  
→ a 4cm egg slice (→ 32 slices)

Dips in 0.20ml H<sub>2</sub>O<sub>2</sub> at 70° overnight & usual Triton -

Results over

scint fl added



# Results

#2-1	3 20cpm	#5-1	184 cpm
2	46	2	0
3	9	3	0
↓	2	↓	↓
79	0	80	0

a.e only at origin in each case